

# Instructions for Use

AMODIA DetectLine

**Basic**

**Molecular test system**  
for the detection of

**Amplification products**

based on the AMODIA-LFD  
(LFD: Lateral-Flow Dipstick)



100

REF:  
ADL-Basic



Distributed by:

AMODIA Bioservice GmbH

Rebenring 31  
D-38106 Braunschweig  
Germany





Tel.: +49 (0) 531-260 17 64  
Fax: +49 (0) 531-260 17 66






E-mail: [info@amodia.de](mailto:info@amodia.de)  
Internet: <http://www.amodia.com>

## Table of Contents

Explanation of Symbols.....	2
Description of the Test.....	3
Materials Supplied, Storage and Stability.....	3
Reagents Required.....	3
Laboratory Instruments and Materials Required.....	3
Warnings and Precautions.....	3
Test Characteristics.....	4
Product Detection.....	4
1.1 Detection.....	4
1.2 Interpretation of results.....	5
Trouble-Shooting.....	6
Method / Test principle.....	6

## Explanation of Symbols

Symbol	Explanation
	Expiry date
	<i>In Vitro</i> Diagnostic Medical Device
	Batch code
REF:	Catalog number
	Storage conditions

Symbol	Explanation
	Consult Instructions for Use
	Consult attended documents
	Package size
	Manufacturer
	Only for evaluation purposes

## Description of the Test

The AMODIA DetectLine Basic is a detection system for double-labeled amplification products.

This test-kit does not contain material to generate those detectable amplification products.

## Materials Supplied, Storage and Stability

Components	Cat-No.	Content (for 25 tests)	Preparation	Storage	Stability / Shelf life
<b>Detektion</b>					
Lateral Flow Dipsticks	LFD01	4 vials a 25 Stk.	ready to use	2 - 8°C*	Until expiration date; 60 days after opening
Chromatographic Buffer	ChB01	2 vials a 10 ml	ready to use	2 - 8°C	Until expiration date; 60 days after opening

\* : LFD vial must be locked tightly ! Storage with opened LFD vial reduces the stability of the LFDs.

### Important note:

Expiration dates should not be exceeded.

## Reagents Required

None

## Laboratory Instruments and Materials Required

- Adjustable pipettes for 10 µl and 200 µl
- Sterile pipette tips with contamination protection (filter tips)
- Single reaction vials or a microwell plate

## Warnings and Precautions

All reagents of this test kit are strictly intended for the specified diagnostic use only. Use by staff, who is especially trained in those methods.

Please adhere strictly to the sequence of processing steps provided by this protocol.

Store all reagents in the original vials at the temperatures indicated on the respective labels. Do not interchange kit components of different lots and assays. Do not use kit components beyond their expiration dates.

Stick to the safety rules for handling kit reagents and sample materials. Especially be aware of the following precautions:

- do not eat, drink or smoke
- wear safety clothes and gloves
- avoid contact with reagents and sample material

Some reagents contain preservation substances against microbial growth, so avoid contact with skin and/or mucous membranes.

Empty vials could be discarded with the normal laboratory waste.

## Test Characteristics

Distributor:	AMODIA Bioservice GmbH
Order-No.:	ADL-Basic
Package size:	100 reactions
Delivery:	from stock in Braunschweig, Germany
Storage:	v. materials supplied, storage and stability

Test time and procedure:


Detection and read-out	approx. 15 minutes
------------------------	--------------------

## Product Detection

### Attention:

- Do not interchange components of different lots.
- Perform the product detection by all means at an area different from the amplification area. (**post-amplification area**).

**Lateral-Flow Dipsticks** (v. fig. 1) are used for the detection of amplification products. They consist (from the bottom) of an area for the chromatographic buffer (white), a spot of sample application (purple), a membrane and an absorption area (white). With the exception of the application area the whole strip is covered by a foil. The strip can be touched at the covered areas. Remarks should be made only on the foil above the absorption area.



**Fig. 1:** Set-up of the Lateral-Flow Dipstick

### 1.1 Detection

	Steps for detection
1.	Take the required number of LFDs from the vial and label them. Only touch the areas covered with foil and use the white area at the end of the LFD for labeling. Close the LFD vial tightly.
2.	For each sample dispense <b>150 µl chromatographic buffer ChB01</b> in single reaction vials or in wells of a (clean) microwell plate.
3.	Pipet <b>5 - 10 µl</b> of the solution containing the double-labeled amplification product onto the application area (purple) at the edge of the foil. A bleeding of the liquid is normal. Incubate for <b>1 min</b> .
4.	Dip the LFDs with the membrane into the chromatographic buffers prepared in step 2 until the application area and the membrane besides the lines are fully discolored (approx. 10-20 minutes). The control line must be visible. Do not read the result before the end of the incubation time. The lines are stable and can be read later.

## 1.2 Interpretation of results


1.	<p><b>Two</b> lines are visible: the test line and the assay control line</p> <p><b>Note:</b> Even a faint test line has to be interpreted as positive. Consider the PCR negative control for a comparison. If needed the test must be repeated for confirmation. Depending on the concentration positive results may be visible even before the entire incubation time is over.</p>	<p>The detection of the double-labeled amplification product is <b>positiv</b>.</p>
2.	<p><b>Only one</b> line gets visible at the position of the assay control line</p> <p><b>Note:</b> The application area should have no violet staining anymore before reading-out the result.</p>	<p>The detection of the double-labeled amplification product is <b>negativ</b>.</p>

The result of a test is only valid if the control line of every sample is stained.

The **positive and negative controls** for the complete test run have to be correct in order to validate the results.

The **positive control** must show a test line **clearly visible**.

The PCR **negative control** must show **no visible test line**. If a stained line occurs, the analysis of **all** the samples tested in parallel must be repeated. If only the signal of the extraction control is positive, a different procedure has to be applied (v. Trouble-Shooting).



**Fig. 2:** Interpretation of the Lateral-Flow Dipstick

The AMODIA DetectLine Basic only gives a qualitative result. The intensity of the stained test line has no direct relation to the number of double-labeled amplification products present in a positive sample.

## Trouble-Shooting

Problem	Possible cause	Recommendation
No assay control line visible	a) Wrong or disfunctional chromatographic buffer ChB01 b) Test strips expired c) Wrong storage of the test strips (moist strips)	Use new chromatographic buffer.  Use new test strips. <b>Store at 2 - 8°C. Close the vials tightly !</b>
All samples and the positive control show a negative signal	a) Amplification not successful b) Labeling of the amplification products not successful	Check the amplification products on a agarose gel. Check the reagents of the labeling.
All samples and controls show a positive signal	a) Contamination of the PCR reaction  b) Chromatographic buffer contaminated with target organism and/or amplification product b) Other reagents contaminated with target organism and/or amplification product	a) Decontaminate the PCR work place and all pipettes. Use new gloves and new pipet-tips from a new box. b) Use fresh chromatographic buffer  c) Use fresh reagents

## Method / Test principle

The detection of the double-labeled amplification products is performed using an immuno-chromatographic assay on a Lateral-Flow Dipstick. These amplification products bind with one label to an antibody which is immobilised on gold particles (gold conjugate). With the diffusion of the chromatographic buffer all gold particles diffuse through the membrane. This membrane is prepared with two lines of different capture molecules. At the first line the capture molecules bind to the second label of the amplification product. So gold particles bound to amplification products are accumulated at this site, thus forming a visible line. Gold particles without amplification product move on. The second line gets visible if the gold conjugate is still intact. This serves as an assay control for the correct performance of the Lateral-Flow Dipstick.

The sensitivity of the AMODIA Lateral-Flow Dipstick is slightly better than a DNA detection with an agarose gel stained with Ethidiumbromide, which is known to detect approx. 10 ng DNA.