

easyStart

All you need to run easyFlow® test-kits

- ▶ all equipment necessary to run easyFlow® test-kits and start right away




- Test-kits:
 - 1x easyFlow® Bacteria (50 rx)
- Devices:
 - 1x incubator
 - 1x mini centrifuge
 - 1x programmable thermal block (pre-installed program)
 - 3x pipettes with variable volumes (20µl, 200µl, 1000µl)
 - 1x rack for liquid culture tubes
 - 1x rack for reaction tubes
 - 1x rack for PCR tubes
- Material:
 - 1x Eugon LT broth liquid culture tubes (100x 9ml)
 - 100x reaction tubes 1,5ml
 - 2x 96 PCR tubes
 - 3x 96 pipette tips 2-200µl
 - 2x 96 pipette tips 50-1000µl


Accelerated Microbial Quality Management

Start with your standard procedure, but be three times faster.

The "easy way":




Standard procedure:



The easy way (~ 25h)

Standard procedure (min.: 72h)



Standard cultivation tests for non-specific microorganisms

*) easyFlow® QC Combo 1 (*S. aureus*, *E. coli*, and *P. aeruginosa*)
easyLAMP® *Candida albicans*

Contact

AMODIA Bioservice GmbH

Rebenring 31
D-38106 Braunschweig
Germany

☎ +49 531 / 260 17 64

📠 +49 531 / 260 17 66

✉ info@amodia.de

🌐 www.amodia.com

