easyFlow[®]

Detect Non-Specified Microorganisms Fast

- test for non-specified microorganisms (ISO 18415) in only one hour of processing
- innovative combination of nucleic acid hybridisation with the lateral-flow technique
- results visible with the naked eye within 60 min
- available for the following species groups:
 - Bacteria
 - Yeast & Mould
 - QC Combo 1 (S. aureus, E. coli, and P. aeruginosa)



60 minutes

Product benefits



- replicating microorganisms only
- nucleic acid based detection
- robust LFD technology




- short hands-on time
- · fast hybridisation protocol



- · no expensive equipment
- · no amplification
- · stable. direct read-out





- · earlier product release
- less consumables




Accelerated Microbial Quality Management

Start with your standard procedure, but be three times faster.







Contact

AMODIA Bioservice GmbH

Rebenring 31 D-38106 Braunschweig Germany

- **4** +49 531 / 260 17 64
- **+**49 531 / 260 17 66
- www.amodia.com

