


easyLAMP[®]

Verify Specified Microorganisms Fast

- ▶ test for specified microorganisms (ISO 18415) in only one hour of processing
- ▶ high sensitivity by isothermal amplification (LAMP) combined with the lateral-flow technique
- ▶ results visible with the naked eye within 60 min
- ▶ available for the following species:
 - *Candida albicans*
 - *Staphylococcus aureus*
 - *Pluralibacter gergoviae*
 - *Pseudomonas aeruginosa*
 - others coming soon



Product benefits



- nucleic acid based detection
- stable, clearly visible detection



- no expensive equipment
- isothermal amplification



- short hands-on time
- fast detection protocol




- earlier product release
- less consumables


Accelerated Microbial Quality Management

Start with your standard procedure, but be three times faster.

The "easy way":




Standard procedure:



The easy way (~ 25h)

Standard procedure (min.: 72 h)



Standard cultivation tests for
non-specific microorganisms

*) easyFlow® QC Combo 1 (*S. aureus*, *E. coli*, and *P. aeruginosa*)
easyLAMP® *Candida albicans*

Contact

AMODIA Bioservice GmbH

Rebenring 31

D-38106 Braunschweig

Germany

☎ +49 531 / 260 17 64

📞 +49 531 / 260 17 66

✉ info@amodia.de

🌐 www.amodia.com

