

DetectLine Basic

- ▶ detection platform for **one** double-labeled molecule

DetectLine Basic *plus*

- ▶ detection platform for **two** double-labeled molecules
- ▶ allows molecular rapid tests with inhibition control

MultiFlow[®]

- ▶ detection platform for up to **four** double-labeled molecules

Product benefits

- standard setup for all kinds of rapid tests:
 - **amplification products** (PCR or isothermal)
 - **labeled antibodies**
- robust detection based on a **lateral-flow strip**
- already includes an **assay control line**
- uses **free labels**
- assay buffer included
- **fast development** (5-10min)



AMODIA

[APPLIED MOLECULAR DIAGNOSTICS]

Develop Your Own Assay

For a molecular assay just use a labeled primer and a labeled probe.
Add an inhibition control to the assay (Basic *plus* and MultiFlow® only).

For an immunological assay just label the antibodies accordingly.

Assay Development Service

If you don't want to develop the assay yourself: We can do that for you.
Just get in contact! Your rapid test may soon be real.

Contact

AMODIA Bioservice GmbH

Rebenring 31

D-38106 Braunschweig

Germany

☎ +49 531 / 260 17 64

📅 +49 531 / 260 17 66

✉ info@amodia.de

🌐 www.amodia.com

